



Corneal Comfort and Protection.

Dissolvable bandage contact lens for use following surgery, injection and traumatic/non-traumatic corneal conditions.

Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times.

 **SOFT SHIELD<sup>®</sup>**  
Collagen Shield

- *Hydrates to 70% water content*
- *Lubricates the ocular surface with collagen while it gradually dissolves*
- *Protects the surface of the eye*
- *Collagen certified suitable for use by EDQM*
- *Comfort of a thin, clear, and pliable collagen barrier*
- *CE approved*

Diameter: 14.5mm Compound base curve: 9mm



*Nourishing and sustaining the eye-care  
community—for the long term.*

Scan Here



Schedule a Knowledge Transfer Session

Call (844) 820-8940

Visit [www.oasismedical.com](http://www.oasismedical.com)

Email [customerservice@oasismedical.com](mailto:customerservice@oasismedical.com)

International [international@oasismedical.com](mailto:international@oasismedical.com)