



# OASIS®

The new SOFT PLUG® Extended Duration 180-T Tapered Canalicular plug is designed with a conical tip for ease of insertion. This medium-term degradable plug lasts up to 6 months in the canaliculus blocking tear drainage through the lacrimal drainage system.

It is an option for patients seeking non-permanent punctal occlusion.




## **SOFT PLUG®** EXTENDED DURATION 180-T Tapered Canalicular Plug

- *Tapered end designed for ease of insertion & placement*
- *Degradable*
- *Lasts up to 180 days*
- *Sterile, insert with forceps*
- *Effective pre and post-surgical use*
- *2mm length*



## The OASIS Punctal Plug Collection



*Nourishing and sustaining the eye-care community—for the long term.*



Scan Here

**View the Punctal Occlusion Research Page**

**Call** (844) 820-8940

**Visit** [www.oasismedical.com](http://www.oasismedical.com)

**Email** [customerservice@oasismedical.com](mailto:customerservice@oasismedical.com)

**International** [international@oasismedical.com](mailto:international@oasismedical.com)

LIT-OV-SPEX180T-S Rev1 04-2023