



Introducing OASIS® Iris Expander.

A surgical solution to gently expand insufficiently dilated pupils during ophthalmic surgery and sustain visibility throughout the procedure.

OASIS® Iris Expander

- *Four pockets that gently cradle the ocular tissue without pinching or clamping*
- *Once removed, the iris returns to its normal shape and function*
- *Also used to manage Floppy Iris Syndrome (IFIS)*
- *Easy removal without potential tearing*
- *Available in sizes 6.25mm & 7.0mm diameter*



Patent No. 8,439,833 B2
Ref # 9625 6.25mm OASIS Iris Expander
Ref #9700 7.00mm OASIS Iris Expander

Nourishing and sustaining the eye-care community—for the long term.

Scan Here



View the Iris Expander Research Page

Call (844) 820-8940

Visit www.oasismedical.com

Email customerservice@oasismedical.com

LIT-OS-IX-S Rev1 04-2023