

OASIS®

The SOFT PLUG Extended Duration 180 Canalicular Plug is a medium-term degradable occluder that lasts up to 6 months, designed to be inserted through the punctal opening into the canaliculus in order to block tear drainage through the lacrimal drainage system. It is an option for patients seeking non-permanent punctal occlusion.




SOFT PLUG®


EXTENDED DURATION 180 Canalicular Plug

- *Degradable*
- *Lasts up to 180 days*
- *Sterile, insert with forceps*
- *Effective pre and post-surgical use*

Available sizes: 0.2mm, 0.3mm, 0.4mm, 0.5mm



The OASIS Punctal Plug Collection



Nourishing and sustaining the eye-care community—for the long term.



Scan Here

View the Punctal Occlusion Research Page

Call (844) 820-8940

Visit www.oasismedical.com

Email customerservice@oasismedical.com

International international@oasismedical.com

LIT-OV-SPEX180-S Rev1 04-2023